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INSTRUCTION MANUAL

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CE 0459

X-Cid®

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CE 0459

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INSTRUCTION MANUAL

CE 0459

This device has been designed and manufactured in accordance with EU Directive 93/42/CE. X-Cid® is an automated device for cleaning handpieces, contra-angles, and turbines.

SAFETY INSTRUCTIONS

- X-Cid® must only be used by qualified personnel to pre-disinfect handpieces, contra-angles, and turbines.
- The instrument must be in good working order. It may not be possible to properly clean a damaged instrument.
- Always wear clean, disposable medical gloves when handling handpieces, contra-angles, turbines and X-Cidol®.
- Only use manufacturer recommended products with the X-Cid®:
 - X-Cidol® 2 detergent; Ref.: 50704005
 - Demineralised water
 - 500I CO2 MICRO-MEGA® oil; Ref.: 50700029
 - The device must be installed by a manufacturer-certified technician following the instructions in the installation manual.
- Do not store the device near a heat source.
- Do not expose the X-Cid® to direct sunlight.
- Do not smoke in close proximity to the device.
- The air exchange rate in the room where the device is installed must be at least 100m³/h.
- In the event the device is not used for four (4) weeks, complete an empty cycle before installing instruments in the treatment tank.
- If detergent comes in contact with the eyes, rinse immediately and thoroughly with water and consult a physician. Consult the X-Cidol® 2 safety data sheet available from your distributor or MICRO-MEGA®.
- Do not disassemble or damage the device.
- Do not place your hands on the hinges when opening the window.
- Make sure the device is solidly seated before use.
- The wall socket must be grounded.
- To unplug the device, remove the power cable.
- Unplug the device and contact your distributor:
 - If the device seems to be operating abnormally (time cycles not respected, abnormal liquid consumption, abnormal noise, etc.)
 - If the liquid is no longer evacuated from the reservoir.
 - If there is a leak at the back of the device.
 - If there is visible damage to the device.
- Contact your distributor or MICRO-MEGA® in order to dispose of the device at the end of its service life in accordance with EU Directive 2002/96/CE.

GENERAL

Introduction

- The X-Cid® automatically pre-disinfects, cleans and lubricates handpieces, contra-angles, and turbines with no disassembly required.
- Treated instruments are clean and in perfect working order, ready to be packaged and sterilised.
- Since the X-Cid® complies with Administrative Circular N° DGS/5C/DHOS/E2/2001/138 of 14 March 2001, patients and health personnel can be assured that there is no risk of cross-infection. In addition, the X-Cid® increases the service life of instruments, reduces maintenance costs, and is compatible with all brands of handpieces, contra-angles, and turbines.

PROCEDURE

■ Pre-disinfection (15 min)

All surfaces in contact with pathogens are treated with X-Cidol® 2, an aldehyde-free detergent/disinfectant solution that complies with current standards. The internal parts of instruments are rotated at low speed to eliminate interfering substances as well as most microorganisms.

■ Rinsing (2 min 30 s)

The rotating instruments are rinsed with demineralised water to evacuate the X-Cidol® 2 used in the previous phase.

■ Cleaning (5 min 30 s)

The instruments are again rotated in X-Cidol® 2. This phase complements the pre-disinfection phase and eliminates all remaining microorganisms, leaving the instruments totally clean.

■ Rinsing (2 min 30 s)

The rotating instruments are rinsed with demineralised water to evacuate the X-Cidol® 2 used in the previous phase.

■ Drying (3 min)

Once clean and rinsed, the instruments are dried using compressed air.


■ Lubricating (2 min 45 s)

The rotating instruments are then lubricated using special oil applied using a mounted spray can. The excess lubricant is removed using compressed air.

DESCRIPTION

- The device is mounted on a stainless steel frame with a self-extinguishing synthetic hull and includes:


1) Treatment Tank:

■ The treatment tank, which has a window, has three stations that can accommodate instrument couplings in any order. Each coupling is held in place by two screws (1). Two types of couplings can be used:

- Type ISO 3964 for handpieces and contra-angles (2).
- Type ISO 9168 (MID) for turbines (3); these couplings can be equipped with quick-connectors for each turbine.

■ The order of the coupling can be changed to accommodate the types of instruments being used.

■ Each station is equipped with four spray nozzles to treat the external surfaces of the instruments (4).

2) Product reservoir:

■ The reservoir is located in the top of the device and is composed of two compartments with independent hinged lids. It has been designed for use with the following products:

- X-Cidol® 2 detergent / disinfect for the pre-disinfection and cleaning phases.
- Demineralised water for the rinsing phases.

■ To fill, the covers of the two compartments are opened and the appropriate quantity of product is added up to the fill line using a funnel, if necessary.

■ Do not exceed the fill line to avoid overflows.

3) Oil spray can:

■ The MM 500 I CO2 oil spray can*, with its cap removed, is placed in its dock on the right hand side of the device. A set of sensors monitors the correct positioning and pressure of the spray can.

▲ Remove cap from spray can.


◀ Place oil spray can in its dock.

* Use only MICRO-MEGA® 500I CO2 oil spray cans (CO2 propellant)

4) Keypad:

■ The control keypad includes the following components:

- A cycle start button
- Three buttons to select the desired stations
- Green, yellow and red LED indicators:
 - Power on (green)
 - Stations selected and cycle underway (yellow)
 - Malfunction (red).

INSTALLATION / USE

■ The X-Cid® must be installed in a well ventilated room. The air exchange rate must be at least 100m³/h per device. The air exchange may be passive or may be assisted by static extractors, air conditioners, or mechanical extractors. The technician will provide the end user with advice in this regard.

■ The X-Cid® is designed to be installed on a horizontal surface, but may be wall mounted using the following method:

- Make sure the wall can support the weight of the device (approx. 8 kg).
- Place two anchors (rawl plugs) for Ø 5 to 5.5 mm screws in the wall 200 mm apart. The screws must be screwed into the anchors so that their heads (Ø 10 maxi) protrude 7 mm from the wall.

■ Empty the reservoirs and remove the 500 l CO2 spray can before moving the device.

■ If the X-Cid® is not used for a year or more, it should be inspected by a certified installer.

EQUIPMENT

■ Equip each station with the following parts, depending on the instruments to be treated:

- ISO 3694 coupling (contra-angles and handpieces).
- ISO 9168 MID coupling (turbines) with suitable quick-connector.

■ Compressed air connection (1)

The device is equipped with a quick-connector for Ø 4x6 flexible tubing. Check that the tubing has a clean cut to ensure an airtight connection.


■ Electrical connection (2)

Plug in the X-Cid® using the power cord provided. Comply with the safety data sheet. The wall socket must be easily accessible and must be grounded.

■ Fluid evacuation (3)

The tubing for evacuating waste fluids (Ø 25 ext.) is directly connected to the waste outlet. All components are biodegradable.

Make sure that the connection is watertight and that the tubing has a sufficient downward slope.

DIRECTIONS FOR USE

- 1) Make sure that the fluids exceed the minimum required levels.
- 2) Turn on the device by pressing the 0/1 button. The green "power on" LED will light up.
- 3) Open the hinged door of the treatment tank and install the instruments to be treated:
 - Insert contra-angles and handpieces completely. You will hear a click. To remove them, press the button and pull.
 - Place turbines on the quick-connect coupling.
- 4) Close the door.
- 5) Press buttons 1, 2, and/or 3 to activate the corresponding treatment stations. The yellow LEDs of the selected stations will blink.

• **Caution ! Instruments should undergo another complete treatment cycle in the event the cycle is voluntarily or involuntarily interrupted in order to ensure an effective treatment.**

- Only select those stations with an attached instrument. Otherwise, cleaning product will be wasted.
- Make sure the stations to be used are correctly selected on the keypad.

6. Press the "start cycle" button.

- During the cycle, the yellow "cycle" and "stations in use" LEDs will light up.
- At the end of a cycle, the "cycle" LED will turn off and the "stations in use" LEDs will blink.
- If the user voluntarily stops the cycle, the "stations in use" LEDs will turn off (which indicates a reset). The user must reselect the stations before starting a new cycle.

Non-lubricating mode

You may also operate the X-Cid® in non-lubricating mode by pressing the "1" button at the same time as the "power on" button.

If the non-lubricating mode is selected, the red LED will light up for 5 seconds when the power is turned on.

MAINTENANCE AND TIPS

The X-Cid® has been designed to clean and lubricate handpieces, contra-angles, and turbines after each patient.

The following tips will help you make the best use of the device and your instruments:

- If your instruments are not new, we recommend a complete internal cleaning to remove accumulated debris prior to a treatment with the X-Cid®.
- Treat your instruments immediately after use so that saliva and blood do not have time to coagulate or dry.
- Once treated, instruments must be packaged for sterilization.
- During long dental procedures (> 15 min) with fast rotating instruments (contra-angles with red rings or turbines), we recommend that you manually lubricate your instruments during the procedure because the air flow tends to eliminate the lubricant.

Daily maintenance

- Before the first daily use of the X-Cid®, verify the fluid levels. In the event you forget this step, the cycle will stop automatically and an alarm will be triggered when fluid levels fall too low.
- Clean the outside of the device with a soft towel dipped in disinfectant. Never use abrasives.

Monthly maintenance

- Select all three stations, operate the device without instrument, and verify that there is a good flow of fluid at each station.

Annual maintenance

- The volumetric flow of the device must be tested one a year by a manufacturer-certified technician.

Fuses

- The fuses are located near the ON/OFF switch. They must be replaced by fuses of the same type (see Characteristics).

CHARACTERISTICS

Definition:

The X-Cid® is an automated pre-sterilization cleaning device.

- Ambient temperature for optimal operation: 10°C to 40°C (50°F to 104°F)
- Ambient humidity for optimal operation: 10% to 80% RH
- Atmospheric pressure for optimal operation: 500 hPa (millibars) to 1,060 hPa (millibars).

Dimensions: H=47cm /L=40cm /W=23cm

Weight: 8 kg

Electrical:

Ref.: 51100003 = 230V +/- 10% 50-60 Hz EU

Ref.: 51100004 = 230V +/- 10% 50-60 Hz CH

Ref.: 51100005 = 230V +/- 10% 50-60 Hz UK

Ref.: 51100007 = 110V +/- 10% 50-60 Hz USA

Ref.: 51100009 = 100V +/- 10% 50-60 Hz Japan

Electric class: I

Power: 38 Watts

Protection :

- Two 0.2 Amp fast blow fuses (5x20 mm) for the 230 V X-Cid®
- Two 0.5 Amp fast blow fuses (5x20 mm) for the 100 and 110 V X-Cid®

Compressed air feed:

- Feed pressure: max.=8 bars / min.=5 bars
- Average flow: 2 l/minute
- Air quality: 1 micron filter; residual water/oil=0.5 ppm at 25°C
- Connection: quick-connect for 4x6 tubing


Evacuation: Tube (Ø 25 mm ext.) connected directly to the waste fluid outlet without siphon.

Treatment fluids:

- Pre-disinfection and cleaning: X-Cidol® 2 detergent; Ref.: 50704005
- Rinsing: Demineralised water
- Lubricating: MICRO-MEGA® 500 I CO2 oil; Ref.: 50700029

INSTRUMENTS TREATED

- The X-Cidol® 2 detergent solution penetrates the components and structures of all commercially available handpieces, contra-angles and turbines, from the bottom to the orifice at the head.


- In case of any doubts, contact your distributor.
- It is essential to carefully examine your old instruments before treating them with the X-Cid®.
- The manufacturer reserves the right to modify the construction and characteristics of the device without prior notice.

WARRANTY

- For any manipulations other than those described above, please contact the service department of your distributor.
- The X-Cid® is guaranteed for **24 months** following the installation date against any manufacturing defects, with the exception of damage caused by improper handling or maintenance, the use of unsuitable products or products that are not recommended by the manufacturer (detergent/disinfectant, demineralised water and oil), or non-compliance with the operating instructions for the device.

DISPLAYS

Indicators

NORMAL OPERATION DISPLAYS

See X-Cid® protocol sheets; Ref. 60300550

• Cycle underway with station 1 selected

• Cycle underway with stations 1 and 2 selected

• Cycle underway with 3 stations selected


• 3 stations selected before cycle start
• End of cycle with 3 stations selected

MALFUNCTION DISPLAYS


See X-Cid® protocol sheets; Ref. 60300550

• Window not properly closed and station 1 selected


• Window not properly closed and stations 1 and 2 selected

• Window not properly closed and all stations selected

• Liquid not present (H₂O or X-Cidol® 2)

• Air not present

• Motor malfunction

• Oil spray can empty or improperly positioned

We recommend that the operator wear a surgical mask and gloves in the event the device malfunctions or the window opens.

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